FV1000 Images Opened With ImageJ

- 1. Download files from OMAL server.
- 2. Open ImageJ.
- 3. In ImageJ: File ---> Open
- 4. Click unzipped file folder to open it.
- 5. Double click to open up image folder you want.
- 6. Click to open up .tif file for that image.
- 7. Repeat step 6 until all .tif files are open
- 7. In menu bar: Image ---> Color ---> Merge Channels.
- 8. Assign proper color to each image. If there is no image for that color channel select *none*. Check Create Composite.
- 9. Click OK.